 SHERA Werkstoff-Technologie GmbH | Espohlstr. 53 | 49448 Lemförde  
Germany | Tel.: +49 (0) 5443 9933-0 | info@shera.de | www.shera.de

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## SHERANAXOS CE0425

Cobalt based dental alloy for prostheses, type 5\*

\*According to ISO 22674:2022

The product is intended for dental laboratories for the creation of products that must be rigid and resistant, single crowns, bridges and superstructures on implants. The alloy SHERANAXOS is distinguished by an outstanding fluidity, which grants to fill even the thinnest details of the frame, down to three tenth of millimeter. Their molecular structure allows to obtain smooth, compact surfaces with little oxide formation

The alloy SHERANAXOS is produced in conformity with standards ISO 22674:2022, it is highly corrosion resistant, it does not contain Nickel and it is completely free of toxic elements like beryllium, indium and gallium.

SHERANAXOS is a product intended to be used by qualified technician.

### Intended purpose

Fabrication of metal plates for RPDs.

### Intended patient population

This alloy is intended to be used for oral rehabilitation in edentulous or partially edentulous patients who suffer from loss of hard dental tissue. There are no limitations related to age, sex, clinical condition or other conditions, other than those possibly determined case by case by the professional user.

### Intended users

The product is intended for professionals of the dental sector, in particular: prosthodontic laboratories, dental technicians, dentists. These professional figures shall have a high degree of specialisation and are perfectly aware of the destination of the product and the correct method of use for the realization of a final product and elaborate different methods to reach the best shape of the final prosthesis. Contact and/or manipulation by other people is not envisaged.

### Instruction for use

#### Modeling

Keep a thickness not less than 0.6 mm, so you can get finished parts up to 0.3mm thick. Avoid elbows, and be careful not to incorporate impurities during waxing.

#### Casting Sprues

Use direct casting sprues of 3.5-4 mm diameter.

#### Casting

Cast SHERANAXOS alloys in pre-heated and cleaned ceramic crucibles used for this alloy only. Each alloy should have its crucible; do not overheat the alloy and do not use the flux.

Reach a temperature between 950°C and 1050°C; maintenance time in pre-heated oven depends on investment characteristics and on the dimension of the cylinder.

Induction casting: start the casting when the cylinders will melt together and just before the superficial crust opens.

Open flame melting: place ingots in the pre-heated ceramic crucible and heat them evenly with circular movements. When ingots have melted, start centrifugal unit.

Use multi-flame welding torches only. Do not use any flux. Indicative values for flame regulation: Acetylene 0,4 bar / Oxygen 2 bar; Propane 0,2 bar / Oxygen 2 bar.

### Manufacturing

Let the cylinder cool down at room temperature, remove the investment and sandblast with aluminium oxide 110 to 150 µm at a pressure of 3-4 bar. Continue manufacturing process using tungsten carbide burs.

### Re-use

The device can be reused until the available material is exhausted.

### Re-use of sprues

The best results are obtained with pure SHERANAXOS. It is nonetheless possible to re-use sprues and cones only once, given that new metal coming from the same lot is added in equal quantity.

### Soldering

We suggest using MAGNUM SALDATURA Co to weld.

### Waste disposal

The processing scrap must be disposed of as special waste in accordance with the EC directives 2008/98/CEE on waste, and 94/62/CEE on packaging and packaging waste and in compliance with national legislation in force on the subject.

### Safety instructions

- Cobalt-Chrome-based alloys can seldom cause dermatitis on sensitive subjects. A Patch-Test is thus advisable.
- Metal dusts and smoke are dangerous for health. Use exhaust fans while casting and polishing.
- The alloy contains Cobalt, element classified as CMR (CAS Number: 7440-48-4).
- Before prosthesis application verify if other metal elements are in patient's oral cavity. Coexistence of different metals can cause a "pile" effect.
- We identify every batch with a number. We recommend to write it down in patient's file to allow its complete traceability.
- This product does not need any special preservation precautions. MESA keeps batch traceability for 15 years.
- The methods of preserving the products made with SHERANAXOS are under the responsibility of professional operators (dental technicians). This product must be stored, transported and used at a temperature between -30°C and +50°C to ensure the integrity of the packaging and the instructions contained therein.
- MESA recommends that the patient be made aware of the possibility for dental alloys to affect MRI results.
- The product is sold non-sterile.

Users shall report any major accident occurred in relation to the use of this medical device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

A summary of safety and clinical performance (SSCP) is available on Eudamed.

Physical and mechanical properties	
Density	8.3 g/cm <sup>3</sup>
Solidus-liquidus temperature	1340-1400°C
Melting point	1450°C
Yield load strength (Rp 0.2)	610 MPa
Percentage elongation at fracture	6%
Modulus of elasticity	200 GPa
Vickers hardness	389 HV10
Colour	White
Thermal expansion coefficient 25-500°C	14·10 <sup>-6</sup> K <sup>-1</sup>
Thermal expansion coefficient 25-600°C	14.4·10 <sup>-6</sup> K <sup>-1</sup>
Ions release after 7 days	0.49 µg/cm <sup>2</sup>
Maximum firing temperature	980°C
Cytotoxicity test according to ISO 10993-5	Passed
Irritation test according to ISO 10993-10	Passed

Symbol	Description	Symbol	Description
	Batch number		Date of manufacture
	Consult instructions for use		
	Manufacturer		CE Marking
	Weight estimate		Non-sterile
	Medical device		Catalogue number
	Temperature limitation		Contains hazardous substances
	Product UDI Code		Distributor

Percentage chemical composition (m/m)			
Co	62	Mo	5
Cr	31		
Others: C, Si, Fe, Mn			