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INSTRUCTIONS FOR USE

- Indication**
Ready-to-use isolating agent for 3D printed resin-models – e.g., against orthodontic & prosthetic resins (cold polymer).
- Processing details: Preparation:**
 - 3D Resins must be fully post cured, cleaned, and dry. (Incomplete polymerized 3D resins may impair the insulating effect.)
 - **Before application of the isolating agent, block out undercuts with wax.** The isolating agent also sticks to the wax.

- Application:**
 - The insulating agent is applied thinly with a clean natural or synthetic hair brush. Applying it several times increases the insulating effect. Allow to dry well between layers. The insulation must no longer stick.
 - Distribute any puddles that may form (e.g. in deep fissures).
 - We recommend a size 8/10 red sable brush.

Advices and tips:

- For large surfaces, e.g. UJ bases, we recommend working in a layering technique: SHERASEPARAT / SHERAiso-3d / SHERASEPARAT / SHERAiso-3d.
- Removal from the model can be facilitated by blocking out e.g. dorsally in areas behind the A-line or the fold-over folds. When removing, a knife can be carefully pressed into the wax and levered.
- For models with large layer thicknesses or from filament printers, spray the blocked-out model with SHERASEPARAT or SHERASEPAL-U in advance and then allow to dry. This perfectly prepares the surface for SHERAiso-3d.
- Keep the bottle well closed after removal.
- To avoid contamination and the introduction of hardened insulation, we recommend placing it in a container. Excess material must not be refilled! (There is a risk that the insulation will gel and become unusable.)
- The blue colouring allows a good visual inspection of the isolated areas and the layer thickness. The dried insulation layer is easy to see.

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