de

en

fr

es

it

pt

da

sv

no

fi

cs

de

en

fr

es

pt

da

sv

no

fi

cs

de

en

fr

es

pt

da

sv

no

fi

cs

de

en

fr

es

pt

da

sv

no

fi

cs

de

en

fr

es

pt

da

sv

no

fi

cs

de

en

fr

es

pt

da

sv

no

fi

cs

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INSTRUCTIONS FOR USE

MODE D'EMPLOI

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BRUKSANVISNING

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Gebruiksaanleitung

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en Instructions for use

Denture base resin, short-term heat cure polymerizate based on methyl methacrylate in powder and liquid form. Free of cadmium, phthalates and BPA.

1. Indications

For the fabrication of full upper and lower dentures using the injection and pressing / tamping techniques. Classification in accordance with DIN EN ISO 20795-1 type 1 class I and in accordance with Regulation (EU) 2017/745, Annex VII, Chapter II: Class Ia for removable and fixed prostheses. Only for intended use by qualified personnel (dental laboratory / dental practice).

2. Contraindication

Patient contact with the unpolymerised acrylic or the individual components is contraindicated. With a confirmed allergy to one or more of the components of the medical device, refrain from using the device. No known interactions.

3. Properties

Simple processing, very good flowability, excellent fit and guaranteed qualitative results, outstanding mechanical values, absolute colour stability, optimal mucosal compatibility thanks to low residual monomer content, tested and certified biocompatibility.

4. Safety instructions

Contains methyl methacrylate; 1,4-butanediol dimethacrylate; dibenzoyl peroxide. Highly flammable liquid and vapour. Causes skin irritation. May cause an allergic skin reaction. May cause respiratory irritation. Toxic to aquatic life with long lasting effects. Avoid breathing dust/fume/gas/mist/vapours/spray. Wear protective gloves/protective clothing/eye protection/face protection. IF ON SKIN: Wash with plenty of water and soap. Observe further safety instructions in the safety data sheet! Finished product: May cause taste irritation or injuries due to breakage of the product. Contact with alcohol, coffee, tobacco smoke or certain medicines may cause discolouration.

5. Processing instructions

Preparation

Clean the cuvette and the teeth carefully by boiling out and then steaming off wax residue (e.g. using SHERAWAX-EX). The plaster model is wetted until no more bubbles rise up so that during pressure polymerisation and when allowing the acrylic to flow or when pressing, no air is forced from the plaster into the acrylic to form bubbles.

SHERADON does not adhere to highly cross-linked artificial teeth. Artificial teeth must be pretreated (roughened) and/or provided with mechanical retention.

Isolation

Plaster surfaces are wetted with two thin layers of alginate isolation (e.g. SHERATEC), ensuring that the first layer is not yet dry when the second is applied. The alginate isolation must be wetted before further processing with acrylic to prevent any white discolourations on the base surfaces.

Mixing ratio 10 : 4

Recommended mixing ratio: 10 g powder with 4 g liquid. Free dosing is possible but a guarantee of the physical material properties is only given if the components are weighed out.

Pressing

Place sufficient liquid in the mixing bowl (made of PE, PP or silicone) and quickly sprinkle over the powder. If free dosing, sprinkle over enough powder to completely bind the liquid in the mixing bowl.

Then thoroughly mix together with a broad spatula (made of PE, PP or metal). Ensure that there are no bubbles. After about 10 minutes* swelling time, tamping time starts (the dough should no longer be sticky). The acrylic dough is then kneaded until homogeneous. The dough must be processed in about a further 30 minutes* (end of the tamping time).

Pressing

Place the acrylic dough in a cooled, handwarm cuvette. Slowly close the cuvette and steadily increase the cuvette pressing force. Maintain the pressing force after closure of the cuvette for 10 min. Sample pressing is not required thanks to the good malleable consistency of the dough.

Polymerisation

Place the cuvette with clamp in boiling water (95°C, about 3 L water for each cuvette) and polymerise for 20 min while boiling. Place very thick items and workpieces with large cross-sections in boiling water, turn off the heat source and wait 15 minutes. Then bring the water to the boil again and polymerise once more for 10 min while boiling. All other polymerisation techniques are of course also possible with this heat cure polymerisate.

Quick facts:

– Mixing ratio: 10 g : 4 g
– Swelling phase: approx. 10 min*
– Working and pressing phase (tamping time): approx 30 min*
– Polymerisation at 95°C: 20 min*

*The times stated refer to a material and room temperature of 23°C. Higher temperatures shorten these times while lower temperatures extend them.

6. Recommendation after completion

Store the prosthesis in water until it is incorporated. As a general rule, 24-hours storage in water is recommended to achieve further breakdown of the residual monomer content (0.8 % as per DIN EN ISO 20795-1) and thus prevent any taste disturbances.

7. Storage

Do not use the material after the expiry date has passed. Do not store above 25°C. Keep out of direct sunlight. Always seal the container well after use. Keep away from heat / sparks / open flames / hot surfaces. Do not smoke.

8. Handling and disposal instructions

To avoid the release of synthetic polymer microparticles, reduce dust emissions as much as possible. In case of excessive dust formation, use extraction with a fine dust filter. Do not allow to enter waste-water systems. Polymer powders can be disposed of with regular household waste. Do not dispose of monomer components in household waste or the sewer system. Dispose of them in accordance with official regulations. Fully polymerized plastic, however, can be disposed of with household waste.

9. Feedback:

Please quote batch number and article number in all correspondence about the product. According to the EU Medical Device Directive, users/patients are obligated to report serious incidents related to a medical device to the manufacturer and the national competent authorities in the country in which they occurred.

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