[de]

[en]

[fr]

[es]

[it]

[pt]

[da]

[sv]

[no]

[fi]

[cs]

[de]

[en]

[fr]

[es]

[it]

[pt]

[da]

[sv]

[no]

[fi]

[cs]

[de]

[en]

[fr]

[es]

[it]

[pt]

[da]

[sv]

[no]

[fi]

[cs]

[de]

[en]

[fr]

[es]

[it]

[pt]

[da]

[sv]

[no]

[fi]

[cs]

[de]

[en]

[fr]

[es]

[it]

[pt]

[da]

[sv]

[no]

[fi]

[cs]

[de]

[en]

[fr]

[es]

[it]

[pt]

[da]

[sv]

[no]

[fi]

[cs]

[de]

[en]

[fr]

[es]

[it]

[pt]

[da]

[sv]

[no]

[fi]

[cs]

[de]

[en]

[fr]

[es]

[it]

[pt]

[da]

[sv]

[no]

[fi]

[cs]

[de]

[en]

[fr]

[es]

[it]

[pt]

[da]

[sv]

[no]

[fi]

[cs]

[de]

[en]

[fr]

[es]

[it]

[pt]

[da]

[sv]

[no]

[fi]

[cs]

[de]

[en]

[fr]

[es]

[it]

[pt]

[da]

[sv]

[no]

[fi]

[cs]

[de]

[en]

[fr]

[es]

[it]

[pt]

[da]

[sv]

[no]

[fi]

[cs]

[de]

[en]

[fr]

[es]

[it]

[pt]

[da]

[sv]

[no]

[fi]

[cs]

[de]

[en]

[fr]

[es]

[it]

[pt]

[da]

[sv]

[no]

[fi]

[cs]

[de]

[en]

[fr]

[es]

[it]

SHERAPRESS-FLEX

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INSTRUCTIONS FOR USE

MODE D'EMPLOI

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BRUKSANVISNING

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[de]

[en]

[fr]

[es]

[it]

[pt]

[da]

[sv]

[no]

[fi]

[cs]

[de]

[en]

[fr]

[es]

[it]

[pt]

[da]

[sv]

[no]

[fi]

[cs]

[de]

[en]

[fr]

[es]

[it]

[pt]

[da]

[sv]

[no]

[fi]

[cs]

[de]

[en]

[fr]

[es]

[it]

[pt]

[da]

[sv]

[no]

[fi]

[cs]

[de]

[en]

[fr]

[es]

[it]

[pt]

[da]

[sv]

[no]

[fi]

[cs]

[de]

[en]

[fr]

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when allowing the acrylic to flow or when pressing, no air is forced from the plaster into the acrylic to form bubbles.

SHERAPRESS-FLEX does not adhere to highly cross-linked artificial teeth. Artificial teeth must be pretreated (roughened) and/or provided with mechanical retention.

Insulation

Plaster surfaces are wetted with two thin layers of alginate insulation (e.g. SHERATEC), ensuring that the first layer is not yet dry when the second is applied.

The alginate insulation must be dry before further processing with acrylic to prevent any white discolourations on the base surfaces.

Mixing ratio 10 : 7

Recommended mixing ratio: 10 g powder with 7 g liquid. Free dosing is possible but a guarantee of the physical material properties is only given if the components are weighed out.

Place the liquid in the mixing bowl (made of PE, PP or silicone) and quickly sprinkle over the powder. If free dosing, sprinkle over enough powder to completely bind the liquid in the mixing bowl. Then thoroughly mix together with a broad spatula (made of PE, PP or metal). Ensure that there are no bubbles.

Processing time

After about 30 seconds* swelling time, the dough can be cast for about 2–3 minutes*. During the casting phase, the acrylic is poured into the matrix. At the beginning of the modelling time, the acrylic becomes firm, no longer flows out of the matrix and can be modelled. The modelling phase is about 20 minutes*.

Polymerisation

After a maximum of 20 minutes, polymerisation takes place in a pressure polymerisation device at a pressure of 2–3 bar and a water bath temperature of 45 to 50°C for approx. 20 minutes.

Quick facts

- Mixing ratio of powder/liquid: 10 g : 7 g
- Swelling phase: approx. 30 s*
- Casting phase: approx. 3 min*
- Plastic modelling phase: approx. 20 min*
- Polymerisation (2–3 bar, 45–50 °C): approx. 20 min**

*The times stated refer to a material and room temperature of 23°C. Higher temperatures shorten these times while lower temperatures extend them.

**The time stated of approx. 20 minutes requires a previous processing time of approx. 20 minutes. The time in the pressure vessel should last up to the 40 minutes after the mixing.

6. Recommendation after completion

Store the prosthesis in water until it is incorporated. As a general rule, 24-hours' storage in water is recommended to achieve further breakdown of the residual monomer content (3,8 % as per DIN EN ISO 20795-1) and thus prevent any taste disturbances.

7. Storage

Do not use the material after the expiry date has passed. Do not store above 25°C. Keep out of direct sunlight. Always seal the container well after use. Keep away from heat / sparks / open flames / hot surfaces. Do not smoke.

8. Handling and disposal instructions

To avoid the release of synthetic polymer microparticles, reduce dust emissions as much as possible. In case of excessive dust formation, use extraction with a fine dust filter. Do not allow to enter wastewater systems. Polymer powders can be disposed of with household waste. Do not dispose of monomer components in household waste or the sewer system. Dispose of them in accordance with official regulations. Fully polymerized plastic, however, can be disposed of with household waste.

9. Feedback:

Please quote batch number and article number in all correspondence about the product. According to the EU Medical Device Directive, users / patients are obligated to report serious incidents related to a medical device to the manufacturer and the national competent authorities in the country in which they occurred.

Fix teeth set up in wax using a plaster or silicone matrix. Then carefully clean the matrix, model framework and the teeth by boiling out any residual wax (e.g. by using SHERAWAX-EX) and then steam off. The plaster matrix is wetted until no more bubbles rise up so that during pressure polymerisation and

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2. **Contraindication**

Patient contact with the unpolymerised acrylic or the individual components is contraindicated. With a confirmed allergy to one or more of the components of the medical device, refrain from using the device. No known interactions.

3. **Properties**

Comfortable range of processing options, very good flowability, optimal mucosal compatibility thanks to low residual monomer content, more homogeneous surface, very low susceptibility to plaque development, easy to finish and polish, absolute colour stability, thanks to catalyst system with no tertiary amines, tested and certified biocompatibility.

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