SHERA-ORTHOMER

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en INSTRUCTIONS FOR USE

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da INSTRUKTIONER TIL BRUG

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no BRUKSANVISNING

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The alginate insulation must be dry before further processing with acrylic to prevent any white discolourations on the base surfaces.

Mixing ratio 10 : 5

Recommended mixing ratio: 10 g powder with 5 g liquid. When working using the dough technique, the powder and liquid can be freely dosed.

Dough technique

Place sufficient liquid in the mixing bowl and quickly sprinkle over the powder. If free dosing, sprinkle over enough powder to completely bind the liquid in the mixing bowl. Then thoroughly mix together with a broad spatula. Ensure that there are no bubbles. After about 30–45 seconds* swelling time, the dough can be worked. The modelling phase is about 10 minutes*.

Spray-on / salt-and-pepper technique

When the alginate insulation is thoroughly dried, the monomer is applied to the plaster model either from a dosing bottle or using a brush and the model is then immediately sprayed with the polymer. Repeat this process until the desired thickness has been achieved. It is particularly important to ensure that only so much liquid is applied that can be absorbed by the powder. Before polymerisation of the plate, the surface must be once again wetted with the liquid and smoothed.

Recommended equipment:

For mixing the components, it is recommended to use monomer-resistant materials – cups, spatulas, spray bottles and nozzles made of PE, PP, silicone or inert metals such as stainless steel.

For the scattering method, a white, type 3 dental plaster is recommended as the model material.

Polymerisation

The polymerisation is carried out in a pressure polymerisation unit with a pressure of 2–3 bar and a water bath temperature of 45°C to 50°C. Polymerisation in a pressure pot, end at the earliest 25 minutes* after starting to mix the powder/liquid together.

Quick facts

- Mixing ratio of powder/liquid: 10 g : 5 g
- Swelling phase: approx. 45 s*
- Plastic modelling phase: approx. 10 min*
- Polymerisation (2–3 bar, 45–50 °C): at least 15 min*

* The times stated refer to a material and room temperature of 23 °C. Higher temperatures shorten these times while lower temperatures extend them.

Recommendation after completion

Store the appliance in water until it is incorporated. As a general rule, 24-hours storage in water is recommended to achieve further breakdown of the residual monomer content and thus prevent any taste disturbances.

- Storage**

Do not use the material after the expiry date has passed. Do not store above 25°C. Keep out of direct sunlight. Always seal the container well after use. Keep away from heat / sparks / open flames / hot surfaces. Do not smoke.

- Handling and disposal instructions**
- To avoid the release of synthetic polymer microparticles, reduce dust emissions as much as possible. In case of excessive dust formation, use extraction with a fine dust filter. Do not allow to enter wastewater systems. Polymer powders can be disposed of with regular household waste. Do not dispose of monomer components in household waste or the sewer system. Dispose of them in accordance with official regulations. Fully polymerized plastic, however, can be disposed of with household waste.

- Feedback:**
- Please quote batch number and article number in all correspondence about the product. According to the EU Medical Device Directive, users / patients are obligated to report serious incidents related to a medical device to the manufacturer and the national competent authorities in the country in which they occurred.

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