

Nachbelichtungsparameter / Post Curing Parameters

SHERAprint Material	Parameter Form Cure V2	Parameter SHERAprint-cure
SHERAprint-model 2.0 color: stone	Preheat OFF 50°C, 3 Min	SHERAprint-model 2.0 stone
SHERAprint-model 2.0 color: smoky white	Preheat OFF 50°C, 3 Min	SHERAprint-model 2.0 smoky white
SHERAprint-model 2.0 color: desert	Preheat OFF 50°C, 3 Min	SHERAprint-model 2.0 desert
SHERAprint-model 2.0 color: aesthetic red	Preheat OFF 50°C, 3 Min	SHERAprint-model 2.0 aesthetic red
SHERAprint-model plaster color: goldenbrown	Preheat OFF 50°C, 3 Min	SHERAprint-model plaster goldenbrown
SHERAprint-model plaster color: grey	Preheat OFF 50°C, 3 Min	SHERAprint-model plaster grey
SHERAprint-cast&press	Preheat OFF 50°C, 3 Min	SHERAprint-cast&press
SHERAprint-splint taff	Preheat OFF Temp OFF, 3 Min	SHERAprint-splint taff
SHERAprint-splint flex	Preheat OFF Temp OFF, 3 Min	SHERAprint-splint flex
SHERAprint-tray clear 2.0	Preheat OFF Temp OFF, 10 Min	SHERAprint-tray clear 2.0
SHERAprint-gingiva 2.0	Preheat OFF 50°C, 3 Min	SHERAprint-gingiva 2.0